



ALMA BEAUTY REMOVE is a highly effective laser hair removal solution

Using an 810 nm "gold standard" diode laser, REMOVE is the **world's best-selling laser hair removal technology**, found in over 5,300 clinics worldwide.

Laser hair removal is a high-end, professional-grade hair removal method and is among the most effective treatments in the industry for achieving smooth and beautiful hair-free skin.

The **combination** of Alma's **810 nm diode laser** with patented **SHR™ technology** makes REMOVE the safest, most reliable and most effective laser hair removal solution available today.



Features

Light Source: Diode Laser 810 nm

Applicators



SPEED™

Half the treatment time

More treatment sessions in less time with a large 2 cm spot size applicator, allowing quick treatment for larger areas, making the treatment fast and most effective.



COMPACT

Easier to Perform

Ergonomic lightweight design makes treatments easier for practitioners and more comfortable for clients.

Specifications:
Spot size 12x10 mm,
up to 10 J/cm²,
up to 10 Hz



FACIAL TIP

Reaching the Hard to Reach

REMOVE's unique facial tip makes hard to reach areas more accessible, including the ears, nostrils and between the eyebrows.

Specifications:
Spot size 6 mm



DISPOSABLE HYGIENIC TIP

Treatment of Intimate Areas

REMOVE includes unique disposable hygienic tips for treating intimate areas.

10 units per package

System

Weight: 40 kg/ 88 lbs

Dimensions: 56 cm/ 22 in (W) x 47 cm/ 19 in (D) x 45 cm/ 18 in (H)

Power Supply: 120/ 230 VAC, 50/ 60 Hz, Single Phase

Benefits

8 years of clinically proven effectiveness

Portable system, for maximum mobility and flexibility

Intuitive, user friendly 15.6" color touch screen

Easy to operate, with a quick learning curve

Dual mode - SHR™ and HR modes

No consumables, no pulse limits or operating costs, yielding a high ROI



Alma Lasers GmbH
Germany

Nordostpark 100-102
90411 Nürnberg, Germany
Tel: + 49 911 / 89 11 29-0
Fax: + 49 911 / 89 11 29-99

info@alma-lasers.de
www.almalasers.com

© 2019 Alma Lasers Inc. All rights reserved.

Alma, Alma Lasers, the Alma Lasers logo, Alma Beauty, REMOVE and SHR In-motion are either registered trademarks or trademarks of Alma Lasers in the United States and/or other countries.

PBBE2705150_02

www.almabeauty.com

| marketing@almabeauty.com

Alma

For You. For Life.



Connect with Alma Beauty